

# USP44 antibody



Catalog Number: 116683

## Product name

USP44 antibody

## Specificity

Human; other species not tested.

## Antibody description

USP44 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse testis tissue. Observed molecular weight by Western-blot: 81kd

## Preparation

This antibody was obtained by immunization of USP44 recombinant protein (Accession Number: NM\_032147). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116683(USP44 antibody) at dilution of 1:500