

# Anti-NUP155 Rabbit antibody



Catalog Number: 167683

## Product name

Anti-NUP155 Rabbit antibody

## Specificity

Human

## Antibody description

NUP155 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human NUP155

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/1000

ICC/IF: 1/200

FC: 1/50

## Validations



Western blot detection of NUP155 in HeLa,A549,HL-60 cell lysates using NUP155 Rabbit mAb(1:1000 diluted).Predicted band size:155kDa.Observed band size:155kDa.