

# HDAC8 antibody

Catalog Number: 111383

## Product name

HDAC8 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HDAC8 Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in HeLa cells, Raji cells. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of HDAC8 recombinant protein (Accession Number: NM\_018486). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111383(HDAC8 antibody) at dilution of 1:500



IP Result of anti-HDAC8 (IP:Catalog No:111383, 4ug; Detection:Catalog No:111383 1:600) with Raji cells lysate 2800ug.