

Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 pAb



Catalog Number: 166383

Product name

Rabbit IgG

Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 pAb

Applications

WB, IHC

Specificity

Human

Dilutions


WB 1:500 - 1:1000

IHC 1:50 - 1:100

Antibody description

Polyclonal antibody to Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220

Validations



Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S246 of human HDAC4

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 pAb

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts from 293 cells untreated or treated with EGF using Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 Antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.

Clonality

Polyclonal

Ig Type