

# Il17a antibody

Catalog Number: 107383

## Product name

Il17a antibody

## Specificity

Human; other species not tested.

## Antibody description

Il17a Mouse Monoclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Il17a recombinant protein (Accession Number: NM-010552). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107383(Il17a Antibody) at dilution of 1:2000