

CLEC3B antibody

Catalog Number: 109383

Product name

CLEC3B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CLEC3B Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Positive IP detected in human plasma tissue. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of CLEC3B recombinant protein (Accession Number: NM_003278). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:109383(CLEC3B Antibody) at dilution of 1:600



IP Result of anti-CLEC3B (IP:Catalog No:109383, 4ug; Detection:Catalog No:109383 1:300) with human plasma tissue lysate 4000ug.