

FSHR antibody

Catalog Number: 110783

Product name

FSHR antibody

Immunogen

[Human FSHR Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

FSHR Rabbit Polyclonal antibody. Positive FC detected in HeLa cells. Positive WB detected in DU 145 cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 71-78 kDa

Preparation

This antibody was obtained by immunization of FSHR recombinant protein (Accession Number: NM_000145). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:110783(FSHR antibody) at dilution of 1:1000



1×10^6 HeLa cells were stained with .2ug FSHR antibody (Catalog No:110783, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.