

IMPAD1 antibody



Catalog Number: 111783

Product name

IMPAD1 antibody

Immunogen

[Human IMPAD1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

IMPAD1 Rabbit Polyclonal antibody. Positive WB detected in 3T3-L1 cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of IMPAD1 recombinant protein (Accession Number: NM_017813). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:111783(IMPAD1 Antibody) at dilution of 1:300