

SYNGAP1 antibody

Catalog Number: 115783

Product name

SYNGAP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SYNGAP1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 148 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006772). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:200-1:2000

IF: 1:10-1:100


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115783(SYNGAP1 antibody) at dilution of



1:500 Immunofluorescent analysis of HepG2 cells, using SYNGAP1 antibody Catalog No:115783 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-SYNGAP1 (IP:Catalog No:115783, 4ug; Detection:Catalog No:115783 1:500) with mouse brain tissue lysate 4000ug.