

# DNAJB2 antibody



Catalog Number: 110093

## Product name

DNAJB2 antibody

## Immunogen

[Human DNAJB2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DNAJB2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A431 cells, HeLa cells, MCF7 cells, mouse brain tissue, mouse heart tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 36-40 kDa

## Preparation

This antibody was obtained by immunization of DNAJB2 recombinant protein (Accession Number: NM\_006736). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110093(DNAJB2 antibody) at dilution of 1:300



IP Result of anti-DNAJB2 (IP:Catalog No:110093, 4ug; Detection:Catalog No:110093 1:1000) with mouse brain tissue lysate 4000ug.