

# NDP Polyclonal Antibody



Catalog Number: 161093

---

## Product name

NDP Polyclonal Antibody

## Specificity

Rat

## Antibody description

Polyclonal antibody to NDP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP\_000257.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions

IHC 1:50 - 1:100

## Validations



Immunohistochemistry - NDP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using NDP antibody at dilution of 1:100 (40x lens).