

Catalog Number: 101093

Product name

Anti-WAP5 / WFDC2 / HE4 antibody

Specificity

Rat WAP5 / WFDC2 / HE4

Antibody description

Rabbit polyclonal to WAP5 / WFDC2 / HE4

Preparation

Produced in rabbits immunized with purified, recombinant Rat WAP5 / WFDC2 / HE4 isoform 2 (Q8CHN3-2; Met 4-Phe 120). WAP5 / WFDC2 / HE4 specific IgG was purified by Rat WAP5 / WFDC2 / HE4 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 µg/mL


This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat WAP5 / WFDC2 / HE4. The detection limit for Rat WAP5 / WFDC2 / HE4 is < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

Validations



WAP5 / WFDC2 / HE4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



WAP5 / WFDC2 / HE4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry