

# HCRTTR1 Polyclonal Antibody



Catalog Number: 166093

## Product name

HCRTTR1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to HCRTTR1

## Preparation

Antigen: A synthetic Peptide of human HCRTTR1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB

## Dilutions

WB 1:500 - 1:1000

## Validations

NCI-H292



Western blot - HCRTTR1 Polyclonal Antibody

Western blot analysis of extracts of NCI-H292 cells, using HCRTTR1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.