

6*His, His-Tag antibody

Catalog Number: 107093

Product name

6*His, His-Tag antibody

Specificity

Tag

Antibody description

6*His, His-Tag Mouse Monoclonal antibody.
Positive IP detected in Transfected HEK-293 cells.
Positive WB detected in recombinant protein.
Positive IF detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000


IP: 1:5000-1:50000

IF: 1:10-1:100

Validations




Western blot of 6*His-tagged fusion protein with anti-6*HIS tag (Catalog No:117294) at various dilutions.



dilutions.

IP Result of anti-6*His, His-Tag (IP:Catalog No:117294, 7ug; Detection:Catalog No:117294 1:10000) with Transfected HEK-293 cells lysate



300ug.

Immunofluorescent analysis of Transfected HEK-293 cells using Catalog No:117294(6*His, His-Tag Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)