

Evi2a antibody

Catalog Number: 110493

Product name

Evi2a antibody

Specificity

mouse; other species not tested.

Antibody description

Evi2a Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of Evi2a recombinant protein (Accession Number: NM_001033711). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110493(Evi2a antibody) at dilution of 1:300