

Human FNDC4/FRCP1 (Fc Tag) recombinant protein



Catalog Number: 502493

General Information

Gene Name Synonym

Fibronectin type III repeat-containing protein 1

Protein Construction

A DNA sequence encoding the human FNDC4 (NP_073734.1) (Met1-Thr167) was expressed with the Fc region of mouse IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 85.4 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 45

Molecular Mass

The recombinant human FNDC4 consists 357 amino acids and predicts a molecular mass of 39.7 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FNDC4 / FRCP1 Protein (Fc Tag) SDS-PAGE