

# BHLHE40 antibody

Catalog Number: 108493

## Product name

BHLHE40 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BHLHE40 Rabbit Polyclonal antibody. Positive IP detected in SGC-7901 cells. Positive WB detected in SGC-7901 cells. Observed molecular weight by Western-blot: 46-50 kDa

## Preparation

This antibody was obtained by immunization of BHLHE40 recombinant protein (Accession Number: NM\_003670). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:117138(BHLHE40 Antibody) at dilution of



1:600

IP Result of anti-BHLHE40 (IP:Catalog No:117138, 4ug; Detection:Catalog No:117138 1:600) with SGC-7901 cells lysate 2200ug.