

JNK2 antibody

Catalog Number: 111893

Product name

Polyclonal

JNK2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

JNK2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of JNK2 recombinant protein (Accession Number: NM_001364608). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111893(JNK2 Antibody) at dilution of 1:300

Store at -20°C. DO NOT ALIQUOT

Clonality