

Anti-ACKR4 antibody



Catalog Number: 175893

Product name

Polyclonal

Anti-ACKR4 antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, IF

Antibody description

Rabbit polyclonal antibody to ACKR4

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ccr11. aa range:291-340

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of various cells using CCRL1 Polyclonal Antibody