

DEPP antibody



Catalog Number: 109893

Product name

DEPP antibody

Specificity

Human; other species not tested.

Antibody description

DEPP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of DEPP recombinant protein (Accession Number: NM_007021). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109893(C10orf10 Antibody) at dilution of 1:600