

NIT2 antibody



Catalog Number: 113193

Product name

NIT2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NIT2 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human brain tissue, human heart tissue. Observed molecular weight by Western-blot: 31kd

Preparation

This antibody was obtained by immunization of NIT2 recombinant protein (Accession Number: NM_020202). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113193(NIT2 antibody) at dilution of 1:500