

BAZ1B Polyclonal Antibody



Catalog Number: 166193

Product name

BAZ1B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BAZ1B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human BAZ1B (NP_115784.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC, IF

Dilutions


WB 1:500 - 1:2000
IHC 1:50 - 1:200
IF 1:50 - 1:200

Validations



Western blot - BAZ1B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BAZ1B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 5s.


Western blot - BAZ1B Polyclonal Antibody

Western blot analysis of extracts from normal (control) and BAZ1B knockout HeLa cells, using BAZ1B antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 1s.

BAZ1B Polyclonal Antibody




Catalog Number: 166193



Immunohistochemistry - BAZ1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using BAZ1B antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BAZ1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse lung using BAZ1B antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - BAZ1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using BAZ1B antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - BAZ1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse stomach using BAZ1B antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - BAZ1B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using BAZ1B antibody at dilution of



1:100 (40x lens).

Immunofluorescence - BAZ1B Polyclonal Antibody

BAZ1B Polyclonal Antibody



Catalog Number: 166193

Immunofluorescence analysis of HeLa cells using

BAZ1B antibody at dilution of 1:100 (40x lens).
Blue: DAPI for nuclear staining.