

Anti-OR4Q3 antibody



Catalog Number: 176193

Product name

Anti-OR4Q3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR4Q3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or4q3. aa range:264-313

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1::; Western Blot analysis of various cells using Olfactory receptor 4Q3 Polyclonal Antibody diluted at 1: 1000