

ERCC1 antibody



Catalog Number: 107193

Product name

ERCC1 antibody

Specificity

Human; other species not tested.

Antibody description

ERCC1 Mouse Monoclonal antibody. Positive WB detected in A431 cells. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of ERCC1 recombinant protein (Accession Number: NM_001369408). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:50-1:500

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:107193(ERCC1 Antibody) at dilution of 1:2000



Immunofluorescent analysis of (10% Formaldehyde) fixed A431 cells using Catalog No:107193(ERCC1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)