

Anti-AMACR (2F3) Mouse antibody



Catalog Number: 168193

Product name

Anti-AMACR (2F3) Mouse antibody

Specificity

Rat

Antibody description

AMACR (2F3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human AMACR(C-terminus) protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/100

Validations



Immunocytochemistry staining of Jurkat cells fixed with -20u2103 Ethanol and using anti-AMACR (C-terminus) mouse mAb (dilution 1:100).



Western blot detection of AMACR(C-terminus) in Rat Liver lysates using AMACR(C-terminus) mouse mAb (1:1000 diluted). Predicted band size: 42KDa. Observed band size: 42KDa.