

# ARIH2 antibody

Catalog Number: 108193

## Product name

ARIH2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARIH2 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HEK-293 cells, HeLa cells, human brain tissue, human colon tissue, human heart tissue, mouse embryo tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 59 kDa

## Preparation

This antibody was obtained by immunization of ARIH2 recombinant protein (Accession Number: XM\_017005532). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108193(ARIH2 antibody) at dilution of 1:500



IP Result of anti-ARIH2 (IP:Catalog No:108193, 4ug; Detection:Catalog No:108193 1:500) with HeLa cells lysate 1200ug.