

# COL4A1 Polyclonal Antibody



Catalog Number: 160593

## Product name

COL4A1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to COL4A1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1445-1669 of human COL4A1 (NP\_001836.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:200 - 1:500

## Validations



Western blot - COL4A1 Polyclonal Antibody

Western blot analysis of extracts of U-87 MG cells, using COL4A1 Antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.