

Anti-NRN1 antibody



Catalog Number: 175593

Product name

Anti-NRN1 antibody

Specificity

Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit

Antibody description

Rabbit polyclonal antibody to NRN1

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human neuritin 51-142/142

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample: Brain(Mouse) lysate at 30ug;;
Primary: Anti-Neuritin at 1:300 dilution;;
Secondary: HRP conjugated Goat-Anti-Rabbit
IgG(bse-0295G-HRP) at 1: 5000 dilution;;
Predicted band size :16 kD; Observed band size :
20 kD