

BRMS1L antibody



Catalog Number: 108593

Product name

BRMS1L antibody

Specificity

Human; other species not tested.

Antibody description

BRMS1L Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, Hela cells. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of BRMS1L recombinant protein (Accession Number: NM_032352). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117240(BRMS1L Antibody) at dilution of 1:600