

KIAA1468 antibody



Catalog Number: 111993

Product name

KIAA1468 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

KIAA1468 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A549 cells. Observed molecular weight by Western-blot: 100-135kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020854). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111993(KIAA1468 antibody) at dilution of 1:500