

SCGBL antibody



Catalog Number: 114993

Product name

SCGBL antibody

Specificity

Human; other species not tested.

Antibody description

SCGBL Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human gliomas tissue, human prostate hyperplasia tissue.

Preparation

This antibody was obtained by immunization of SCGBL recombinant protein (Accession Number: NM_001025591). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:114993(SCGBL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114993(SCGBL antibody) at dilution of 1:50 (under 40x lens)