

# PCSK2 Polyclonal Antibody



Catalog Number: 165993

## Product name

PCSK2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PCSK2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human PCSK2 (NP\_002585.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



## Western blot - PCSK2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PCSK2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.