

Anti-ADGRG5 antibody

Catalog Number: 175993

Product name

Anti-ADGRG5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ADGRG5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gpr114. aa range:91-140

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of COLO205 cells using GPR114 Polyclonal Antibody