

Ascrinvacumab reference antibody



Catalog Number: 176993

Product name

Ascrinvacumab reference antibody

Specificity

Human ACVRL1 / ALK-1

Antibody description

Anti-ACVRL1 / ALK-1 Antibody (Ascrinvacumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-ACVRL1 / ALK-1 Antibody (ascrinvacumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-ACVRL1 / ALK-1 Antibody (ascrinvacumab) is more than 95%, determined by SEC-HPLC.