

Hepsin antibody



Catalog Number: 111293

Product name

Hepsin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Hepsin Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001110252). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111293(Hpn antibody) at dilution of 1:400