

# NOL3 antibody

Catalog Number: 113293

## Product name

NOL3 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

NOL3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IF detected in HepG2 cells. Positive WB detected in mouse skeletal muscle tissue, HeLa cells. Observed molecular weight by Western-blot: 30-33kd




## Preparation

This antibody was obtained by immunization of NOL3 recombinant protein (Accession Number: NM\_001276312). Purification method: Antigen affinity purified.

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113293(NOL3 antibody) at dilution of 1:800

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113293(NOL3 antibody) at dilution of 1:100 (under 40x lens)

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunofluorescent analysis of HepG2 cells, using

# NOL3 antibody



Catalog Number: 113293

---

NOL3 antibody Catalog No:113293 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit  
IgG (red).