

# B4GALT4 Polyclonal Antibody



Catalog Number: 164293

## Product name

B4GALT4 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to B4GALT4

## Preparation

Antigen: A synthetic peptide of human B4GALT4

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - B4GALT4 Polyclonal Antibody

Western blot analysis of extracts of mouse skeletal muscle, using B4GALT4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.