

# SLAIN2 antibody



Catalog Number: 115293

## Product name

SLAIN2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SLAIN2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of SLAIN2 recombinant protein (Accession Number: NM\_020846). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115293(SLAIN2 Antibody) at dilution of 1:300