

# Rabbit anti GST-Tag pAb



Catalog Number: 166293

## Product name

Rabbit anti GST-Tag pAb

## Antibody description

Polyclonal antibody to Rabbit anti GST-Tag

## Preparation

Antigen: Recombinant protein of human GST-Tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:1000 - 1:2000

IF 1:20 - 1:50


IP 1:20 - 1:50

## Validations



Western blot - Rabbit anti GST-Tag pAb

Western blot analysis of GST-tagged protein using GST-tag antibody at different dilution. Each lane was loaded with 5 ug protein. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Blocking buffer: 3% nonfat dry milk in TBST.



TBST.

Western blot - Rabbit anti GST-Tag pAb

Western blot analysis of over-expressed GST fusion protein using GST-tag antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.