

FOXC2 antibody



Catalog Number: 110693

Product name

FOXC2 antibody

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

FOXC2 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, A431 cells, rat spleen tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of FOXC2 recombinant protein (Accession Number: NM_005251). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110693(FOXC2 antibody) at dilution of 1:1000