

PCDHA2 antibody

Catalog Number: 113693

Product name

PCDHA2 antibody

Immunogen

[Human PCDHA2 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PCDHA2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, human heart tissue, Raji cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 80 kDa

Preparation

This antibody was obtained by immunization of PCDHA2 recombinant protein (Accession Number: NM_031496). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113693(PCDHA2 antibody) at dilution of 1:800



IP Result of anti-PCDHA2 (IP:Catalog No:113693, 3ug; Detection:Catalog No:113693 1:300) with HeLa cells lysate 4650ug.