

STAT5B antibody



Catalog Number: 115693

Product name

STAT5B antibody

Specificity

human; other species not tested.

Antibody description

STAT5B Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 91kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_012448). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115693(STAT5B antibody) at dilution of 1:1500