

# Anti-FGFR2 antibody



Catalog Number: 176693

## Product name

Anti-FGFR2 antibody

## Dilutions

WB: 1:2,000

## Specificity

ICC: 1:200

Human

## Validations

## Antibody description

Rabbit polyclonal antibody to FGFR2

## Preparation

This antigen of this antibody was recombinant protein within human fgfr2/cd332 aa 600-750 / 821.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC, IHC-P, FC




Fig1: Western blot analysis on cell lysates using anti- CD332 rabbit polyclonal antibodies.




Fig2: ICC staining CD332 in SK-BR-3 cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS and counterstained with DAPI in order to highlight the nucleus (blue).