

# THY1 Polyclonal Antibody



Catalog Number: 163393

## Product name

THY1 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to THY1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 31-81 of human THY1 (NP\_006279.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - THY1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using THY1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - THY1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using THY1 antibody at dilution of 1:100 (40x lens).