

Anti-ZNF695 antibody



Catalog Number: 176393

Product name

Anti-ZNF695 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZNF695

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human znf695.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293T cells using ZNF695 Polyclonal Antibody diluted at 1: 2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).