

CLIP-115/CYLN2 antibody



Catalog Number: 109393

Product name

CLIP-115/CYLN2 antibody

Clonality

Polyclonal

Immunogen

[Human CLIP2 Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

CLIP-115/CYLN2 Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 110-120 kDa


Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of CLIP-115/CYLN2 recombinant protein (Accession Number: NM_032421). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109393(CLIP2 Antibody) at dilution of 1:3000