

# Mouse XEDAR/EDA2R (Fc Tag) recombinant protein



Catalog Number: 500793

## General Information

### Gene Name Synonym

X-linked ectodysplasin-A2 receptor

### Protein Construction

A DNA sequence encoding the mouse EDA2R (Q8BX35) (Met1-Thr138) was expressed, fused with the Fc region of human IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Met

### Molecular Mass

The recombinant mouse EDA2R /Fc is a disulfide-linked homodimer. The reduced monomer

comprises 379 amino acids and has a predicted molecular mass of 42.4 kDa. The apparent molecular mass of the protein is approximately 49 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse XEDAR / EDA2R Protein (Fc Tag) SDS-PAGE