

# MTP18 antibody

Catalog Number: 112793

## Product name

MTP18 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MTP18 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 16-18 kDa

## Preparation

This antibody was obtained by immunization of MTP18 recombinant protein (Accession Number: NM\_016498). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112793(MTP18 antibody) at dilution of 1:300



IP Result of anti-MTP18 (IP:Catalog No:112793, 4ug; Detection:Catalog No:112793 1:500) with mouse heart tissue lysate 4000ug.