

Syntaxin 16 antibody



Catalog Number: 115793

Product name

Syntaxin 16 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Syntaxin 16 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, Jurkat cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of Syntaxin 16 recombinant protein (Accession Number: NM_003763). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:115793(STX16 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115793(STX16 antibody) at dilution of 1:50 (under 10x lens)