

Anti-eIF3B Rabbit antibody



Catalog Number: 166793

Product name

Anti-eIF3B Rabbit antibody

Specificity

Human

Antibody description

eIF3B Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human eIF3B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/50

Validations



Western blot detection of eIF3B in K562, HeLa cell lysates using eIF3B Rabbit pAb(1:1000 diluted). Predicted band size: 93kDa. Observed band size: 117kDa.